

Single-use Minimally Invasive Surgical Products

10mm Single-use Clip Applier



5mm Single-use Clip Applier




Trocar



Irrigation/Suction Kit



Economic I/S set



Auto-Retrieval Series



Single-use Clip Applier

SOLIDCLIP™

10mm CA16 Series

Multi-fire Titanium Clip Applier
CE, TFDA Certified (FDA to be obtained)



Anti-drop Design

of the
Jaw Structure.

**Pitch-Loading
Clip feeding Design**

contributes a
Lighter & Stable
Clips Ligating
Experience.

Interlock Design

of the trigger
when
running out of
the clips.

Patented Pitch-Loading design

Product	Model No.	Description	Dimensions
Clip Applier	CA-165	Single-use, L Size, 10 Clips	10x330mm
Clip Applier	CA-166	Single-use, L Size, 15 Clips	10x330mm
Clip Applier	CA-167	Single-use, M/L Size, 15 Clips	10x330mm
Clip Applier	CA-168	Single-use, M/L Size, 20 Clips	10x330mm
Clip Applier	CA-169	Single-use, M/L Size, 10 Clips	10x330mm



PMDA, CE MDD, TFDA Certified
(FDA, CE MDR to be obtained soon)

Single-use Clip Applier

SOLiDCLIP™

5mm CA56 Series

Multi-fire Titanium Clip Applier



Anti-drop Design

of the
Jaw Structure.

**Pitch-Loading
Clip feeding Design**

contributes a
Lighter & Stable
Clips Ligating
Experience.

Interlock Design

of the trigger
when
running out of
the clips.

Patented Pitch-Loading design

Model No.	Description	Dimensions
CA-567	Single-use, M/L Size, 10 Clips	5x330mm
CA-568	Single-use, M/L Size, 15 Clips	5x330mm

PMDA, TFDA Certified
(FDA, CE to be obtained soon)

Trocar

- **Bladeless Dilating Trocar**

- Exclusive funnel providing optimal target for smooth instrument exchange.
- Tapered threads which provide smooth insertion and superior hold in abdominal wall.
- Easily removed funnel cap for specimen removal & rapid de-sufflation.

Product	Model No.	Description	Dimensions
Trocar	AST-PB80P0	Bladeless dilating trocar & rigid cannula with stopcock Set	5x80mm
Trocar	AST-PB10K0		5x100mm
Trocar	AST-MD10A0		10x100mm
Trocar	AST-MD70A0		10x70mm
Trocar	AST-ME10A0		12x100mm
Trocar	AST-ME70A0		12x70mm



Universal seal build-in for widespread instrument exchange.



Thread design for better stability




Unique blade shield for smoother insertion

Cadence Suction / Irrigation set

- Supports Pistol and Trumpet Grips
- Comfort, controlled suction and flow are essential ingredients to a surgical procedure
- Reliable suction performance & Intuitive design
- Precision irrigation & effortless flow

Product	Model No.	Description	Dimensions
Cadence SI	18-001L	3.7m tubing, dual spike	5x330mm
Cadence SI	18-002L	3.7m tubing, dual spike	5x450mm
Cadence SI	18-003L	3.7m tubing, single spike	5x330mm
Cadence SI	18-004L	3.7m tubing, single spike	5x450mm



The Multigate Cadence® Suction Irrigator delivers effortless management of the operating field through precision and flowing ergonomic design.

FEATURES Precision and flow

Comfort, control and flow are essential ingredients to a surgical procedure. Multigate's Cadence® Suction Irrigator delivers the effortless flow, ergonomic comfort, controlled suction and precision you need to perform at your best.

- ✓ Reliable suction performance
- ✓ Precision irrigation
- ✓ Supports pistol and trumpet grips
- ✓ Leak free
- ✓ Surgeon comfort
- ✓ Intuitive design

Leak free: Simplified internal mechanism eliminates leaks.

Reduced fatigue: Ergonomic hand-piece for firm but comfortable instrument grasp reduces surgeon fatigue.


Intuitive: Icons, colours and grooved suction button ensure suction and irrigation functions are easy to identify.

Precision: Ø5mm precision tip provides constant, precise irrigation flow without compromising suction.

Configurable: Dual shaft attachment options allow for quick and secure change between pistol and trumpet grip.



Received for outstanding design and innovation. Find out more: Multigate.com.au/GDA



Economic Suction / Irrigation set

CE
2460



- Quick probe exchange for 10mm probe, diathermy probe & extra length demand
- Phthalate free of PVC tubing set
- Ergonomically handle designed

Model No.	Description	Dimensions
ASI-A20000	premier suction irrigation set with dual spike	5x330mm
ASI-A25000		5x450mm
ASI-A27000		5x330mm



Exchangeable probe



Color coded trumpet valves



Low torque tubing & DEHP free

Best MIS Choice

Insufflation Tubing set

- **99.99% filtration efficiency**

- Filter designed to withstand high gas flow rates and work with high-flow insufflator
- 4 type of connector & compatible with various insufflator connectors.



Low torque tubing & DEHP free



High efficient Filter



Model No.	Description
AIF-4000	Insufflation tubing set with luer lock adaptor
AIF-4300	Insufflation tubing set with push fitting adaptor
AIF-4400	Insufflation tubing set with O-ring adaptor
AIF-4600	Insufflation tubing set with 15/22 adaptor
AIF-6000	Insufflation tubing set with luer lock adaptor
AIF-6300	Insufflation tubing set with push fitting adaptor
AIF-6400	Insufflation tubing set with O-ring adaptor
AIF-6600	Insufflation tubing set with 15/22 adaptor
AIF-6700	Insufflation tubing set with Y connector



luer connector



Push fit



O-ring connector



15/22 connector



Y connector



Auto-Retrival Series

CE
2460



- **Lap Bag with Auto Retrieval Pouch**

- Pouch Self-opening designed and keep round shape of pouch
- U shape pouch and available for 4 sizes
- V shape pouch and available for 2 sizes
- It can be detached pouch from spring-load frame and suture tie for closing bag

Model No.	Description	Dimensions
ALB-CUD0	U shape detachable - Lap Bag with auto retrieval	4"x4"(200ml)
ALB-CUE0	U shape detachable - Lap Bag with auto retrieval	5"x5"(300ml)
ALB-CUF0	U shape detachable - Lap Bag with auto retrieval	6"x6"(800ml)
ALB-CUG0	U shape detachable - Lap Bag with auto retrieval	7"x7"(1350ml)
ALB-AV30	V shape detachable - Lab Bag with auto retrieval	3"x6"(170ml)
ALB-AV50	V shape detachable - Lab Bag with auto retrieval	5"x7"(260ml)



Ergonomic handle

Strong Nylon strip for closed pouch

Best MIS Choice

Lap Bag

- **Model Number:**
 - ALB-GV30
- **Description:**
 - Lap bag with 10mm introducer, 3"x6"(175ml) elastic pouch with suture tie.
- **Quantities:** 3000~5000 pcs
- **Unit Prices:** USD\$6.0 (EXW)
- **Delivery Time:**
 - 1 month after receiving order & payment (valid till March. 1st, 2023)

