



# SOLIDCLIP™ 5 mm

## Single-use Clip Applier



Patented Pitch-Loading design



Anti-drop Design

of the  
Jaw Structure.



Pitch-Loading  
Clip feeding Design

Lighter & Stable  
Ligating Force.



Interlock Design

when running out  
of the clips.

Model					Preloaded Clip Number	Color Indicator			Dimensions
	Clip Size	W	mm H	L		Clip Number			
CA-568	M/L	3.5	8.2	9.0	15 Clips				5 x 330 mm
						1st ~ 12th	13th ~ 15th	No Clips	

**Medscope Biotech Co., Ltd.**

2F, No. 8, Keyi Street, Miaoli County, 350 Taiwan, R.O.C.

Tel: +886-37-580055 www.medscope.com.tw, Email: globalsales@medscope.com.tw

# SOLID CLIP™



## Medscope Laparoscopy



### **Medscope Biotech Co., Ltd.**

2F, No. 8, Keyi Street, GuangYuan Technology Park, Zhunan Town, 35059 Miaoli County, Taiwan

Tel: +886 (037) 580055 Fax: +886 (037) 585050

E-mail: [globalsales@medscope.com.tw](mailto:globalsales@medscope.com.tw)

[www.medscope.com.tw](http://www.medscope.com.tw)